## User Manual

Document No  
*Doküman No*  
TF03-IFU03-02

Issue Date  
*Yayın Tarihi*  
03.06.2022

Revision No  
*Revizyon No*  
01

Revision Date  
*Revizyon Tarihi*  
11.08.2022

### CORTEX EXTRACTION CANNULA *CORTEX EXTRACTION KANÜL*

Company Name <i>Şirket Adı</i>	:	MERAN TIP ÜRÜNLERİ TEKNOLOJİ SAN.TİC.LTD.ŞTİ
Company Representative <i>Firma temsilcisi</i>	:	Cemil BAYTAŞ / General Manager
Address <i>Adres</i>	:	Ferhatpaşa Mah. Arif Nihat Sok. A Blok No:4/1A ÇATALCA/ İSTANBUL/TÜRKİYE
Phone <i>Telefon</i>	:	+90 216 557 8 755
Web <i>Web</i>	:	<a href="http://www.merantip.com">www.merantip.com</a>

#### **Product Structure and General Features:**

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The Cortex Extraction Cannula is designed for the removal of the lens fragments in the cortex by cataract surgeries.. Cannula tip is made of 304 Stainless Steel. The connection area at the back of the cannula is made of polycarbonate raw material. The products are available for single use and sterile. The cannulas to be used may differ depending on the operation area and the surgeon's choice. Models are available in various diameters and sizes in order to be suitable for every eye structure.

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#### **Product Models:**

##### **Ürün Modelleri:**

Model and size information of the Cortex Extraction Cannula products manufactured by us are shown in **TF03.TB.01 Annex-1 Product Models Table** with pictures.

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#### **Intended Use:**

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Cortex Extraction Cannulas are disposable and used for the removal of the lens fragments in the cortex by cataract surgeries.

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**Revision Date**  
*Revizyon Tarihi*  
11.08.2022

### **Indications:**

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Cortex Extraction Cannulas are sterile conducive equipments, which are used to inject or reject fluids for removal of lens cortex in ophthalmic surgical operations.

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### **Patient Population:**

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Patients with various eye conditions, especially cataracts.

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### **Contraindications:**

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It should not be used other than ophthalmic surgery.

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It should not be used on infected eyes.

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### **Complications:**

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If the tip of the Cannula is not sharp, it may damage the tissue of the eye.

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After operation mild chemosis and redness can occur.

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Intraoperative serous occurrence.

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### **Side Effects:**

#### **Yan Etkiler:**

There are no known side effects of the product.

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### **Product Users:**

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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







- During the operation, BSS (Balanced Salt Solution) fluid is injected with the Cortex Extraction Cannula (which is suitable according to the physician's preference).

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- The capsular bag is gently removed via suction to prepare the area for cortex removal.

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
**Revision No**  
*Revizyon No*  
01

**Revision Date**  
*Revizyon Tarihi*  
11.08.2022

- Ophthalmic Viscoelastic Device (OVD) is injected onto the lens to protect the tissue against any damage. The cortex is aspirated using the Cortex Removal Cannula. Thus, the process is completed for the following operation steps.

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### **Factors That Can Affect the Success of The Operation:**

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- The user's experience and experience in the operation,  
*Kullanıcının operasyondaki deneyimi*
- Product selection in suitable type and size for the patient to be applied,  
*Uygulanacak hastaya uygun tip ve ebatta ürün seçimi,*
- Eye structure of the patient,  
*Hastanın göz yapısı,*
- Failure of the patient to have the necessary controls after the operation.  
*Failure of the patient to have the necessary controls after the operation.*

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Revision No  
*Revizyon No*  
01

Revision Date  
*Revizyon Tarihi*  
11.08.2022

### **Product Storage and Shipping Conditions**

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Protect the product against direct sunlight during shipment and storage.

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Keep the product dry in any environment.

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Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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**Revision No**  
*Revizyon No*  
01

**Revision Date**  
*Revizyon Tarihi*  
11.08.2022

### **Medical Equipment Used with Products:**

#### ***Medical Equipment Used with Products:***

The products are used with injectors during surgery.

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### **Accessory:**

#### ***Aksesuar:***

The product does not have any accessories.

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### **Shelf Life:**

#### ***Raf Ömrü:***

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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### **Packaging, Labeling and Handling:**

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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The packaging sample of our products is as seen in the photos below.

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