User Manual

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IRRIGATION AND AIR INJECTION CANNULA *IRRIGATION AND AIR INJECTION KANÜL*

Company Name <i>Şirket Adı</i>	: MERAN TIP ÜRÜNLERİ TEKNOLOJİ SAN.TİC.LTD.ŞTİ
Company Representative <i>Firma temsilcisi</i>	: Cemil BAYTAŞ / General Manager
Address <i>Adres</i>	: Ferhatpaşa Mah. Arif Nihat Sok. A Blok No:4/1A ÇATALCA/ İSTANBUL/TÜRKİYE
Phone <i>Telefon</i>	: +90 216 557 8 755
Web <i>Web</i>	: www.merantip.com

Product Structure and General Features:

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The Irrigation and Air Injection Cannula is designed for the use in cataract surgeries and their purpose is to inject various liquids and air on the eye and anterior chamber. Cannula tip is made of 304 Stainless Steel. The connection area at the back of the cannula is made of polycarbonate raw material. The products are available for single use and sterile. The cannulas to be used may differ depending on the operation area and the surgeon's choice. Models are available in various diameters and sizes in order to be suitable for every eye structure.

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Product Models:

Ürün Modeli:

Model and size information of the Irrigation and Air Injection Cannula products manufactured by us are shown in **TF03.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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Irrigation and Air Injection cannulas are disposable and usually used in cataract surgeries and their purpose is to inject various liquids and air on the eye and anterior chamber. It may also be used for removal of air, fluids & viscoelastic substances during an ophthalmic surgery.

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Indications:

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Irrigation and Air Injection cannulas are used to inject air, fluids and viscoelastic substances into eyes during ophthalmic surgeries.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used other than ophthalmic surgery.

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It should not be used on infected eyes.

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Complications:

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If the tip of the Cannula is not sharp, it may damage the tissue of the eye.

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After operation mild chemosis and redness can occur.

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Intraoperative serous occurrence.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training.

Users must have operational experience.

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Instruction for Use:

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Irrigation and Air Injection Cannulas are restricted for use by surgeons trained in ophthalmic surgery.

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➤ The Cannulas should be opened from the marked are without disruption of the sterility.

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➤ Cannula should be connected to the Equipment by a sterile person according to usage of the product.

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➤ Depending on the use, air, fluids or viscoelastic substances are injected into eye or removed from eye.

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➤ Dispose of used products in accordance with established hospital protocol for hazardous waste.

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Factors That Can Affect the Success of The Operation:

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

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Keep the product dry in any environment.

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Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

Medical Equipment Used with Products:

The products are used with injectors during surgery.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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
Revision No
Revizyon No
01

Revision Date
Revizyon Tarihi
02.06.2023

Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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