User Manual

Manuale d'uso

Document No	Issue Date	Revision No	Revision Date
<i>No documento</i>	<i>Data di rilascio</i>	<i>No revisione</i>	<i>Data revisione</i>
TF03-IFU03-12	03.06.2022	02	02.06.2023

VITREORETINAL SURGERY CANNULA CANNULA PER CHIRURGIA VITREORETINICA

Company Name	:	MERAN TIP ÜRÜNLERİ TEKNOLOJİ SAN.TİC.LTD.ŞTİ
<i>Nome azienda</i>		
Company Representative	:	Cemil BAYTAŞ / General Manager
<i>Rappresentante dell'azienda</i>		
Address	:	Ferhatpaşa Mah. Arif Nihat Sok. A Blok No:4/1A
<i>Indirizzo</i>	:	ÇATALCA/ İSTANBUL/TÜRKİYE
Phone	:	+90 216 557 8 755
<i>Telefono</i>		
Web	:	www.merantip.com
<i>Web</i>		

Product Structure and General Features:

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The Vitreoretinal Surgery Cannula is designed for the use in very different place and purpose by vitreoretinal surgical procedures. Cannula tip is made of 304 Stainless Steel. The connection area at the back of the cannula is made of polycarbonate raw material. The products are available for single use and sterile. The cannulas to be used may differ depending on the operation area and the surgeon's choice. Models are available in various diameters and sizes in order to be suitable for every eye structure.

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Product Models:

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Model and size information of the Vitreoretinal Surgery Cannula products manufactured by us are shown in **TF03.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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Used in very different place and purpose by vitreoretinal surgical procedures. These are;

- Viscoelastic fluid administration across the posterior intraocular lens surface
- Injection of vital dyes and perfluorocarbon liquids during small-gauge vitrectomy surgery
- Performing vitrectomy in endophthalmitis and trauma.

Manuale d'uso

Document No	Issue Date	Revision No	Revision Date
<i>No documento</i>	<i>Data di rilascio</i>	<i>No revisione</i>	<i>Data revisione</i>
TF03-IFU03-12	03.06.2022	02	02.06.2023

- Clear up residual peripheral vitreous

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Indications:

Indicazioni:

Vitreoretinal Surgery Cannulas are used in vitreoretinal surgeries in multiple purposes.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used other than ophthalmic surgery.

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It should not be used on infected eyes.

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Complications:

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If the tip of the Cannula is not sharp, it may damage the tissue of the eye.

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After operation mild chemosis and redness can occur.

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Intraoperative serous occurrence.

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Side Effects:

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There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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Instruction for Use:

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Vitreoretinal Surgery Cannulas are restricted for use by surgeons trained in ophthalmic surgery.

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Factors That Can Affect the Success of The Operation:

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The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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Single Use / Reusability Status:

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The products are for single use only. There is absolutely no reusability.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

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Keep the product dry in any environment.

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Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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The products are used with injectors during surgery.

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Accessory:

Accessori:

The product does not have any accessories.

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User Manual

Manuale d'uso

Document No	Issue Date	Revision No	Revision Date
No documento	Data di rilascio	No revisione	Data revisione
TF03-IFU03-12	03.06.2022	02	02.06.2023

Shelf Life:

Validità:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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
Revision No
No revisione
02

Revision Date
Data revisione
02.06.2023

The packaging sample of our products is as seen in the photos below.

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Data revisione

02.06.2023

Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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