# User Manual

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Document No  
*Doküman No*  
TF08-IFU08-03

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Revision No  
*Revizyon No*  
03

Revision Date  
*Revizyon Tarihi*  
01.06.2023

**MUSCLE HOOK**  
***KAS KANCASI***

**Company Name** : MERAN TIP ÜRÜNLERİ TEKNOLOJİ SAN. TİC. LTD. ŞTİ.  
*Firma Adı*  
**Authorized Representative** : Cemil BAYTAŞ / General Manager  
*Firma Yetkilisi*  
**Address** : Ferhatpaşa Mah Arif Nihat Sok. A Blok No:4/1AÇatalca/İstanbul /  
*Adres* : TÜRKİYE  
**Phone** : +90 216 557 87 55  
*Telefon*  
**Web** : [www.merantip.com](http://www.merantip.com)  
*Web*

### **Product Models:** **Ürün Modelleri:**

Model and size information of the Muscle Hook products produced by us are shown in **TF08.TB.01 Annex-1 Product Models Table** with pictures.

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There are three model of the product. They are “Muscle Hook”, “Muscle Hook(Round)” and “Muscle Depressor(T)”.

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### **Muscle Hook;** **Kas Kancası:**

It is used to move the eye muscles up and down, right and left during the operation in eye surgery. In addition, Muscle Hooks in Eye Surgery are designed to elevate the muscle tissue for dissection and re-attachment in performing strabismus or other procedures with ocular muscles.

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### Muscle Hook(Round);

#### *Kas Kancası(Yuvarlak);*

It is used to move the eyeball to the desired area during eye surgery or to ensure that the eyeball does not move during the operation.

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### Muscle Depressor(T);

#### *Kas Depresörü(T)*

Muscle Depressor is used for controlling eye position. The product is used by the operator to depress the sclera of a patient's eye for either rotation or immobilization of the globe of the eye during examination.

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### Indications:

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Muscle Hooks are generally used to move the eye muscles. In addition, it is used to provide or restrict the movement of the eyeball during the operation.

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### Patient Population:

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Patients with various eye conditions, especially cataracts.

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### Contraindications:

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It should not be used outside of ophthalmic surgery.

It should not be used on infected eyes.

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Revision Date  
*Revizyon Tarihi*  
01.06.2023

### **Complications:**

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After operation mild chemosis and redness can occur.

Intraoperative serous occurrence.

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### **Side Effects:**

#### *Yan Etkiler:*

There are no known side effects of the product.

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### **Product Users:**

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The product should be used by specialist physicians who have completed the necessary medical training.

Users must have operational experience.

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### **Usage Method:**

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- It is used by the surgeon to control the movement of the eye muscles during the operation. After the necessary incisions are made on the eye surface, it is inserted under the muscle layer. It is directed by the surgeon to control the movements of the eye until the end of the operation.

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

Revision No  
**Revizyon No**  
03

Revision Date  
**Revizyon Tarihi**  
01.06.2023

### **Factors That Can Affect the Success of the Operation:**

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### **Shelf Life:**

#### Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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### **Product Storage and Shipping Conditions**

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Protect the product against direct sunlight during shipment and storage.

Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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### **Medical Equipment Used with Products:**

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The products are used with no additional medical equipment.

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### **Accessory:**

#### **Aksesuar:**

The product does not have any accessories.

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### **Packaging, Labeling and Handling:**

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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


Revision No  
**Revizyon No**  
03

Revision Date  
**Revizyon Tarihi**  
01.06.2023

The packaging sample of our products is as seen in the photos below.

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Revision No  
**Revizyon No**  
03

Revision Date  
**Revizyon Tarihi**  
01.06.2023

### Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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