# User Manual

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Document No  
*Doküman No*  
TF09-IFU09-01

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Revision No  
*Revizyon No*  
00

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### Ophthalmic Gases *Oftalmik Gazlar*

Company Name : MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ  
Authorized Representative : Cemil BAYTAŞ / General Manager  
Address : İstanbul Trakya Serbest Bölgesi, Ferhatpaşa SB Mah. Ali Rıza Efendi Cad. No: 22/201 Çatalca/İSTANBUL/TÜRKİYE  
Phone : +0212 786 6253  
Web : [www.merantip.com](http://www.merantip.com)

#### Product Models:

##### Ürün Modelleri:

MMT BRAND <i>MMT MARKA</i>	MERAN BRAND <i>MERAN MARKA</i>	ÜRÜN ADI <i>PRODUCT NAME</i>
<i>SF6-10</i>	<i>M10 1000</i>	Sf6 10 cc Sulfur hexafluoride
<i>SF6-30</i>	<i>M10 1010</i>	Sf6 30 cc Sulfur hexafluoride
<i>SF6-50</i>	<i>M10 1020</i>	Sf6 50 cc Sulfur hexafluoride
<i>C3F8-10</i>	<i>M10 2000</i>	C3f8 10 cc Octafluoropropane
<i>C3F8-30</i>	<i>M10 2010</i>	C3f8 30 cc Octafluoropropane
<i>C3F8-50</i>	<i>M10 2020</i>	C3f8 50 cc Octafluoropropane

#### Intended Use:

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Ophthalmic gases are categorized as Sulphur Hexafluoride(SF6) and Perfluoropropane (Octafluoropropane(C3F8)) separately. A sterile plastic graduated syringes mounted with a sterilizing 0.22 µm filter is connected to aluminum can where ophthalmic gas is kept. Ophthalmic gases are used in pars plana vitrectomy procedures where a retina hole or tear is repaired and act to provide a long-term tamponade, or plug, of a retinal hole/tear and allows re-attachment of the retina to occur over the following several days post-op.

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
### **Indications:**

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Indications of ophthalmic gases are given below:

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Revision No  
*Revizyon No*  
00

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### **Patient Population:**

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There are no restrictions on the patient population. The choice of appropriate patient is the responsibility of the surgeon. However, the products are generally used in patients who will undergo retinal detachment surgery.

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### **Contraindications:**

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- Nitrogen protoxide anaesthesia should not be used when ophthalmic gases are injected.

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- Ophthalmic gases should not be used in the areas other than ophthalmology.

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### **Complications:**

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- In large retinal tears, small gas bubbles may pass into the subretinal space and expand the detachment.

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- During gas injection with topical anesthesia, the patient should be informed immediately if the image becomes gray and disappears, and arterial pulsation must be followed during the injection.

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- Central retinal artery occlusion may occur.

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- Potential side effects would be increase of intraocular pressure, glaucoma, lens opacities, bullous keratopathies, subretinal gas, zonular dehiscence or phacodonesis, new and enlarged tears, dislocation of IOL and must be considered by the surgeon with regards to the expected benefits.

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### **Side Effects:**

#### **Yan Etkiler:**

The complications described in the literature are those linked to vitreo-retinal surgery and particularly:

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- It is known that it can cause Proliferative Vitreoretinopathy and new tear formation by causing cortical vitreous compression on the retinal surface.

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- If people with large perfluorocarbon gas bubbles in the eye experience sudden drops in pressure ((elevation, travel by plane, diving with or without isobar caisson, hyperbar caisson treatment), it causes the gas bubble to expand.

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### **Product Users:**

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The products should be used by experienced surgeons who are experienced in the surgical technique.

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### **Surgical Technique:**

#### Cerrahi Teknik:

*There are a few issues that need to be considered in order to obtain and use the gas as pure as possible. The first of these is the evacuation of the air of the syringe. For this purpose, pure gas can be taken after the injector is filled once or twice and emptied out. After the gas is taken, it should be used immediately. Otherwise, gas may leak out of the plastic injector by diffusion during standby.*

*It is extremely important to loosen all tractions before the liquid air exchange to be performed after vitrectomy. Otherwise, air will escape under the retina during the air exchange and the retina will not be soothed. The air under the retina is still evidence of traction. In this case, the fluid should be switched to immediately and the retina should be mobilized by removing the tractions. Full liquid air exchange is important for the formation of gas bubbles of sufficient size. For this reason, the fluids accumulating in the posterior pole and under the retina should be aspirated repeatedly for at least 10 minutes from the beginning of the air exchange. If gas exchange will follow air exchange, gas injection is started in a controlled manner after one of the sclerotomy sites is fully closed and the other suture is passed, and when the air-gas mixture is used, the air in the eye is considered to be completely replaced with gas.*

*Ensure that the central and peripheral vitreous is carefully removed before injecting the gas. During gas injection, monitoring of correct pupillary irrigation and possible optic nerve compression is essential; In these cases, the injection should be discontinued immediately.*

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### **INSTRUCTIONS FOR USE:**


#### ***KULLANIM TALİMATI***

It is necessary to prepare the mixture with air, before proceeding with the infusion procedure; this is done directly in the syringe, by proceeding as follows:

- connect the filter to the syringe from the female luer inlet fitting side
- fit the filter on the aluminum can
- push the aluminum can to fill the syringe with a small amount of gas (approx. 5ml)
- empty the syringe and refill it with enough gas to cross the line marked as 20%
- push the plunger to position (1) at the suggested mixing value (20%), to expel any excess gas
- suck in the air, bringing the plunger to position (2) indicated by 100%
- never remove the filter from the syringe during these operations
- remove filter once the syringe is full and proceed to use.

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### **Sterility/Sterilization Method:**

#### ***Sterilite/Sterilizasyon Metodu:***

The products are sterilized with Ethylene Oxide gas and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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### **Shelf Life:**

#### **Raf Ömrü:**

Shelf life of the product is determined as 4 years as a result of stability study.

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### **Duration of the Product in the Body:**

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Ophthalmic gases may remain in eyes for 3 weeks. There is no need to perform an additional surgery to remove ophthalmic gases.

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### **Product Storage and Shipping Conditions:**

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Protect the product against direct sunlight during shipment and storage.

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Keep the product dry in any environment.

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Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Document No  
Doküman No  
TF09-IFU09-01

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Revision No  
Revizyon No  
00

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### **Medical Equipment Used with Products:**

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Ophthalmic Gases are composed of the following components;

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Following the Instructions for Use, luer-lock cannula should be attached to syringe for the injection of ophthalmic gases.

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### **Accessory:**

#### ***Aksesuar:***

The product does not have any accessories.

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### **Packaging, Labeling and Handling:**

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Packaging materials are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It is sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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