HYDRODISSECTION CANNULA
HYDRODISSECTION KANÜL

Company Name Şirket Adı	:	MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve Fatma Baytaş
Company Representative Firma temsilcisi	:	Cemil BAYTAŞ / General Manager
Address Adres	:	İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE
Phone Telefon	:	+90 212 786 6253
Web Web	:	www.merantip.com

Product Structure and General Features:

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The Hydrodissection Cannula is designed for the use in the cataract surgery to effectively separate layers of the nucleus and cortex. Cannula tip is made of 304 Stainless Steel. The connection area at the back of the cannula is made of polycarbonate raw material. The products are available for single use and sterile. The cannulas to be used may differ depending on the operation area and the surgeon's choice. Models are available in various diameters and sizes in order to be suitable for every eye structure.

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Product Models:

Ürün Modeli:

Model and size information of the Hydrodissection Cannula products manufactured by us are shown in **TF03.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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Hydrodissection Cannulas are disposable and used in the cataract surgery to effectively separate layers of the nucleus and cortex. Hydrodissection Cannulas also used for aspiration of the nucleus out of the capsular bag and distribution of fluids around the nucleus.

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Indications:

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Hydrodissection Cannulas are used to separate layers of nucleus and cortex to aspirate nucleus in cataract surgeries.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used other than ophthalmic surgery.

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It should not be used on infected eyes.

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Complications:

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If the tip of the Cannula is not sharp, it may damage the tissue of the eye.

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After operation mild chemosis and redness can occur.

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Intraoperative serous occurrence.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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Instruction for Use:

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Hydrodissection Cannulas are restricted for use by surgeons trained in ophthalmic surgery.

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Factors That Can Affect the Success of The Operation:

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Sterility/Sterilization Method:

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The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

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Keep the product dry in any environment.

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Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

Medical Equipment Used with Products:

The products are used with injectors during surgery.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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
Packaging, Labeling and Handling:

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

Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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