LASIK CANNULA
LASIK KANÜL

| | | |
|--|---|--|
| Company Name Şirket Adı | : | MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve Fatma Baytaş |
| Company Representative Firma temsilcisi | : | Cemil BAYTAŞ / General Manager |
| Address Adres | : | İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE |
| Phone Telefon | : | +90 212 786 6253 |
| Web Web | : | www.merantip.com |

Product Structure and General Features:

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The Lasik Cannula is designed for keeping the operation surface wet thanks to its perforated structure and it is used to position the lasik flap in lasik surgeries. Cannula tip is made of 304 Stainless Steel. The connection area at the back of the cannula is made of polycarbonate raw material. The products are available for single use and sterile. The cannulas to be used may differ depending on the operation area and the surgeon's choice. Models are available in various diameters and sizes in order to be suitable for every eye structure.

The Lasik Cannula is designed for keeping the operation surface wet thanks to its perforated structure and it is used to position the lasik flap in lasik surgeries. Cannula tip is made of 304 Stainless Steel. The connection area at the back of the cannula is made of polycarbonate raw material. The products are available for single use and sterile. The cannulas to be used may differ depending on the operation area and the surgeon's choice. Models are available in various diameters and sizes in order to be suitable for every eye structure.

Product Models:

Ürün Modeli:

Model and size information of the Lasik Cannula products manufactured by us are shown in **TF03.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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In Lasik surgeries, Lasik Cannulas keep the operation surface wet thanks to its perforated structure and it is used to position the lasik flap. Multiple ports provide multi-directional irrigation.

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Indications:

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Lasik Cannulas are used to keep operation area wet and to position lasik flap.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used other than ophthalmic surgery.

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It should not be used on infected eyes.

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Complications:

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If the tip of the Cannula is not sharp, it may damage the tissue of the eye.

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After operation mild chemosis and redness can occur.

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Intraoperative serous occurrence.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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Instruction for Use:

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Lasik Cannulas are restricted for use by surgeons trained in ophthalmic surgery.

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Factors That Can Affect the Success of The Operation:

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

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Keep the product dry in any environment.

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Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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
Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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The packaging sample of our products is as seen in the photos below.

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Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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