DCR and Lacrimal Intubation Set
DCR ve Lakrimal Entübasyon Seti

Company Name <i>Firma Adı</i>	:	MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve Fatma Baytaş
Authorized Representative <i>Firma Yetkilisi</i>	:	Cemil BAYTAŞ / General Manager
Address <i>Adres</i>	:	İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE
Phone <i>Telefon</i>	:	0212 786 6253
Web <i>Web</i>	:	www.merantip.com

Product Models:

Ürün Modelleri:

Model and size information of the DCR and Lacrimal Intubation Set products produced by us are shown in **TF01.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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
DCR and Lacrimal Intubation Sets are used for liquid transfer and to unlock the lacrimal canal between eye and nose. Products are sterile and single use conductive equipments. Sample models and explanations for DCR and Lacrimal Set are listed below.

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Lacrimal Intubation Set:

Lakrimal Entübasyon Seti:

Lacrimal Bicanalicular Intubation Set contains two flexible stainless steel probes (0,40-0,90mm diameters) with an olive shaped tip and the other ends joined by a length of silicone tubing. Various sizes are available.

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Curved Probe

Kavisli Prob

DCR Set contains two flexible stainless steel probes (0,30-0,90mm diameters) with an olive shaped tip and the other ends joined by a length of silicone tubing. Probes come in two types: Straight and Curved. Various sizes are available. DCR sets are specially designed for Dacryocystorhinostomy surgeries.

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Indications:

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DCR and Lacrimal Intubation Sets are used when;

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-Where probing of the lacrimal system has not cured the tearing

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-Where there are definite obstructions of the lacrimal system that would close if they were not kept open

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-Repair of injuries involving cut canaliculi

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-During a dacryocystorhinostomy (in particular DCR sets are used).

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Patient Population:

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There is no restriction on patient population. The appropriate patient selection is the responsibility of the surgeon. The products are generally used in patients with obstructed tear ducts.

Contraindications:

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DCR and Lacrimal Intubation is contraindicated:

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Complications:

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As with any surgery, there are risks associated with the procedure used or the development of the initial pathology. Possible complications associated with implantation of DCR and Lacrimal Intubation Sets include, but are not limited to:

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Possible Complications During Surgery:

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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Operational Technique:

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- The patient is brought to the operation area and prepared for surgery. Upper punctum is dilated with Punctal Dilator and Probe device.

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

- After the punctum canal is expanded, the dilator is removed from the operation site. After that, the appropriate product in the DCR and Lacrimal Intubation Set is selected by the surgeon. After the selection, the product is placed in the hole previously created.

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- One of the probes is passed through the upper punctum and then pulled down through the nasal cavity with the help of Retrieval Hook.

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- Then, the length of the silicon tube in the patient's nasal cavity and its tightness between the patient's upper and lower punctum's are checked manually. After the tightness and length of the silicone tube are adjusted, the silicone tubes are tied each other.

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- Finally, the silicon tubes are fixed to the patient's nasal walls and the operation is completed by cutting the long parts of the silicon tubes.

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Factors That Can Affect the Success Of The Operation:

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Sterility/Sterilization Method:

Sterilite/Sterilizasyon Metodu:

The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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During the operation, the punctum of patients with small lacrimal punctum is expanded with the help of a punctal dilator.

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Accessory:

Aksesuar:

Some sets are offered to the user with the retrieval hook accessory.

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
Packaging, Labeling and Handling:

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


Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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

The packaging and boxing of the Lacrimal Intubation Set is as follows.

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Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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