Document No
Doküman No
TF04-MMT-IFU04

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

INTRAVITREAL INJECTION KIT
İNTRAVİTRAL ENJEKSİYON KİTİ

Company Name
Firma Adı : MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve Fatma Baytaş
Authorized Representative
Firma Yetkilisi : Cemil BAYTAŞ / General Manager
Address
Adres : İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE
Phone
Telefon : 0212 786 6253
Web
Web : www.merantip.com

Product Models:

Ürün Modelleri:

Model and size information of the Intravitreal Injection Kit products produced by us are shown in **TF04.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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Intravitreal Injection Kit products are used during Ophthalmic Surgery, to assist before and during intravitreal injection and to reduce the risk of infection.

The content of kit, intended use and usage area is detailed below.

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It is a self-holding ophthalmic surgical instrument used to open the eyelids during the ophthalmological examination or procedure. The metal speculum produced enables the surgeon to work safely in a wide angle by stretching the eyelids of the patient and stopping involuntary blinking.

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➤ Drape with Pouch (40X40);
Kılıflı Örtü (40X40);

It is used to prevent contamination that may occur in the patient during ophthalmic surgery (cataract, intraocular, refractive procedures etc.) and to safely collect the fluids formed in the operation area. Liquids are collected in the Pouch on the product.

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➤ Table Drape (70X80);
Masa Örtüsü (70X80);

Surgical table drapes are used to prevent contact with the patient and surfaces during surgery and to protect peripheral surfaces, equipment and patient from contamination.

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➤ Cotton Swabs;
Pamuklu Cubuklar;

During ophthalmic operations such as cataract surgery, it is used to absorb and remove the discharge from the eye in order for the surgeon to perform easy operation. This product is made of cotton raw material.

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➤ Gauze;
Gazlı Bez;

Medical gauze is placed in the operation area after ophthalmic surgery and used to protect the eye until the wound heals.

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➤ 20G and 30 G Cannula
20G ve 30G Kanül

Ophthalmic Cannulas are designed according to the intended place of use and convenient to the surgical operation and is a sterile surgical conducive equipment, which is using for aspiration, irrigation and tissue contact.

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Injector is used for injection of various fluids during ophthalmic surgery.

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It is used to identify and mark the area where the fluid to be given into the eye during ophthalmic surgery will be given.

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It is worn by operating room staff during ophthalmic surgery and is used to protect the surgical wound from contamination.

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Indications:

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Indications for intravitreal injections may include: diabetic retinopathy; macular degeneration; ocular infection; retinal swelling; retinal vascular disease, and choroidal neovascular membrane.

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Patient Population:

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There are no restrictions in the patient population. Proper patient selection is the surgeon's responsibility. However, the products are used in operations such as diabetic retinopathy, macular degeneration, ocular infection, retinal swelling, retinal vascular disease, choroidal neovascular membrane.

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Contraindications:

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It should not be used outside of ophthalmic surgery.

It should not be used on infected eyes.

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It should not be used on infected eyes.

Complications:

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Intraocular Pressure may occur.

Endophthalmitis (intraocular infection) may occur.

Cataracts may occur.

Retinal detachment may occur if the injection is made through the retina.

If the injection site is too close to the limbus, the needle may damage the lens.

Subconjunctival hemorrhage may occur.

Vitreous hemorrhage may occur.

Retinal artery occlusion may occur.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training.

Users must have operational experience.

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Surgical Technique:

Cerrahi Teknik:

- Intravitreal Injection Kit is brought to the area where the operation will be performed and all products are revealed.

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- Topical anesthesia is performed to perform the operation. 10% povidone iodine is applied to the periorcular area, eyelid, eyelashes and the areas where the sterile cover will come into contact.

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- The Drape with Pouch is placed on the patient's face. The perforated part is adjusted to be on the patient's eye.

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





Factors That Can Affect the Success of The Operation:

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions:

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Protect the product against direct sunlight during shipment and storage.

Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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There is no medical equipment used with our products.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Packaging, Labeling and Handling:

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The packaging sample of our products is as seen in the photos below.

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