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Product Models:

Ürün Modeli

Model and size information of the Laser Drape Film products produced by us are shown in **TF07.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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Laser Drape Film is used in ophthalmic laser operations to protect hygiene, to ensure safety against infection risk for patients and controls fluids out of patients' body.

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Indications:

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Laser Drape Film is used in surgical operations in order to prevent patients from infection which may arise from personnel and operation environment.

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Patient Population:

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Patients who undergo ophthalmic laser operations.

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Contraindications:

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Laser Drape Film should not be used in operations such as general surgery, dental surgery, cardiovascular surgery other than ophthalmic laser surgery.

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Complications:

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If the patient is allergic to the product material, there will be irritation or redness in the areas where the product is in contact with the patient's tissues.

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Side Effects:

Yan Etki

There are no known side effects of the product.

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Product Users:

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The product should be used by medical personnel who attend ophthalmic laser operations.

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Usage Information:

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➤ Before using the drape set, check the package to make sure there are no torn holes.

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➤ Take the product into the sterile area and open it from top or bottom of the bag.

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➤ Unfold and pull off the label which is located around the round hole(s).

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➤ Adhere round hole(s) of laser drape film around patient's eye considering the arrow(s) which indicates direction(s) of pouch(es) on the label located on the round hole by aligning the hole to patient's eye.

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➤ If one pouch exists, put the pouch right downward for right eye and left downward for left eye. If two pouches exist, put pouches right and left side of patient's face downward.

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➤ At the end of the surgical procedure, collect the drapes, place them in a medical waste bag and take them to the medical waste area for disposal.

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Sterility/Sterilization Method:

Sterilite/Sterilizasyon Methodu

The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Shelf Life:

Raf Ömrü:

Shelf life of the product is determined as 4 years as a result of stability study.

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Product Storage and Shipping Conditions

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The ambient temperature where the products are stored must be in the range of +5 / +35 ° C and the humidity at maximum 70% rh.

The products should be stored in a way not to be exposed to direct sunlight and high above the ground to prevent contact with water.

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Medical Equipment Used with Products:

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The products are used with no additional medical equipment.

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Accessory:

Aksesuar

The product does not have any accessories.

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Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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
Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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