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SURGICAL SEAT ARMREST COVER
Cerrahi Çift Koltuk Kılıfı

Company Name <i>Firma Adı</i>	:	MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve Fatma Baytaş
Company Representative <i>Firma Yetkilisi</i>	:	Cemil BAYTAŞ / General Manager
Address <i>Adres</i>	:	İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE
Phone <i>Telefon</i>	:	0212 786 6253
Web <i>Web</i>	:	www.merantip.com

Product Models:

Ürün Modelleri:

Model and size information of the Surgical Seat Armrest Cover products produced by us are shown in **TF07.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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Surgical Seat Armrest Cover is used in all surgical operations to prevent infections, to protect hygiene, to ensure safety against infection risks which may arise from surgical seat armrests.

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Indications:

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Surgical Seat Armrest Cover is used in surgical operations in order to prevent patients from infection and contamination which may arise from personnel and operation environment.

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Patient Population:

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Patients who undergo surgical operations.

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Contraindications:

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Apart from preventing infection, our products are used for radiation etc. It should not be used to protect against such substances.

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Complications:

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If the patient is allergic to the product material, there will be irritation or redness in the areas where the product is in contact with the patient's tissues.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by medical personnel who attend surgical operations.

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Sterility/Sterilization Method:

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The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Shelf Life:

Raf Ömrü:

Shelf life of the product is determined as 4 years as a result of stability study.

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Product Storage and Shipping Conditions

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The ambient temperature where the products are stored must be in the range of +5 / +35 ° C and the humidity at maximum 70% rh.

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The products should be stored in a way not to be exposed to direct sunlight and high above the ground to prevent contact with water.

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Medical Equipment Used with Products:

Medical Equipment Used with Products:

The products are used with no additional medical equipment.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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