User Manual

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Product Models:

Ürün Modelleri:

Model and size information of the Muscle Hook products produced by us are shown in **TF08.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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There is one model of the product. It is "Punctal Dilator and Probe (KHAN)".

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Punctal Dilator and Probe (KHAN):

Punctal Dilatör ve Prob (KHAN)

It is used with the DCR product to open or widen the tear ducts during ophthalmic surgery. Thanks to its bi-directional structure, it provides easy use with one hand for the user.

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Indications:

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Punctal Dilator and Probes are generally used with the DCR product to open or widen the tear ducts during ophthalmic surgery. Thanks to its bi-directional structure, it provides easy use with one hand for the user.

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Patient Population:

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It is used in patients with congestion in tear ducts.

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Contraindications:

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It should not be used outside of ophthalmic surgery. It should not be used on infected eyes.

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Complications:

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After operation mild chemosis and redness can occur. Intraoperative serous occurrence.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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Usage Method:

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- The sharp part of the product is placed in the eye canal of the patient by the surgeon. Then, while the product is pushed into the channel, it is turned around and the channel widening is made.

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
Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

- After this procedure, fluid is introduced into the opened canal with the lacrimal cannula by the surgeon.

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Factors That Can Affect the Success of the Operation:

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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The products are used with no additional medical equipment.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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
Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

The packaging sample of our products is as seen in the photos below.

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Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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