Document No
Doküman No
TF08-MMT-IFU08-06

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

OPHTHALMIC EYE SHIELD
OFTALMİK GÖZ KALKANI

Company Name <i>Firma Adı</i>	:	MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve Fatma Baytaş
Authorized Representative <i>Firma Yetkilisi</i>	:	Cemil BAYTAŞ / General Manager
Address <i>Adres</i>	:	İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE
Phone <i>Telefon</i>	:	0212 786 6253
Web <i>Web</i>	:	www.merantip.com

Product Models:

Ürün Modelleri:

Model and size information of the Ophthalmic Eye Shield products produced by us are shown in **TF08.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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There is one model of the product. It is “Universal Eye Shield”.

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Universal Eye Shield;

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It is placed on the eye socket to protect the eye after ophthalmic surgery. Thanks to its perforated structure, it provides air around the eye.

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Indications:

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Ophthalmic Eye Shields are generally used for the protect of the eye after ophthalmic surgery.

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Revision Date
Revizyon Tarihi
20.11.2023

Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used outside of ophthalmic surgery.

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Complications:

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Cross contamination can occur when the same product is used in more than one patient.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training.

Nurses or health personnels can use the products.

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Factors That Can Affect the Success of the Operation:

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Sterility/Sterilization Method:

Sterilite/Sterilizasyon Metodu:

The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions

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
Protect the product against direct sunlight during shipment and storage.Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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There is no medical equipment to be used with products.

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Accessory:**Aksesuar:**

The product does not have any accessories.

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

The packaging sample of our products is as seen in the photos below.

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