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Product Models:

Ürün Modelleri:

Model and size information of the Pediatric Rop Kit products produced by us are shown in **TF08.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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There is one model of the product. It is "Pediatric Rop Kit".

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Pediatric Rop Kit:

Pediyatrik Rop Kiti:

This kit is a combination of Muscle Hook and Speculum products. In this direction, the purpose of use of the product covers the functions of two products. The intended use of the kit is, move the eye muscles during the operation and keep open the eyelids during the ophthalmological examination applied in pediatric patients.

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Indications:

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Pediatric Rop Kits are generally used to keep the eyelid open and to move the eye muscles during the operation.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used outside of ophthalmic surgery.

It should not be used on infected eyes.

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Complications:

After operation mild chemosis and redness can occur.

Intraoperative serous occurrence.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Revision Date
Revizyon Tarihi
20.11.2023

Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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
Usage Method:

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- This kit is a combination of Muscle Hook and Speculum products. So the usage method of the product is a combination of Muscle Hook and Speculum products too.

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
- It is used by the surgeon to control the movement of the eye muscles during the operation. After the necessary incisions are made on the eye surface, it is inserted under the muscle layer. It is directed by the surgeon to control the movements of the eye until the end of the operation.

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During the operation, the selection of the appropriate size and shape of the product should be made by the surgeon.

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Factors That Can Affect the Success of the Operation:

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Sterility/Sterilization Method:

Sterilite/Sterilizasyon Metodu:

The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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There is no medical equipment to be used with products.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in

User Manual

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Document No
Doküman No
TF08-MMT-IFU08-08

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
Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

desired dimensions. Leak-proof barrier allows the product to remain sterile.

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
Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

The packaging sample of our products is as seen in the photos below.

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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