User Manual

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Product Models:

Ürün Modelleri:

Model and size information of the Silicone Injection-Extraction Set products produced by us are shown in **TF08.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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There are two model of the product. They are “Silicone Injection Kit” and “Silicone Injection-Extraction Set”.

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Silicone Injection Kit:

Silikon Enjeksiyon Kiti:

It is used to supply pressurized air to the silicon injector. The kit is connected to the injector. It is used as a connecting element to deliver the air it takes from other medical devices to the injector.

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Silicone Injection-Extraction Set;

Silikon Enjeksiyon-Ekstraksiyon Seti:

It is used to deliver the silicone used during ophthalmic surgery to the operation area or to remove it from the operation area. This set does not connect to the injector.

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Indications:

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Silicone Injection-Extraction Sets are generally used for deliver the silicone used during ophthalmic sur-gery to the operation area or to remove it from the operation area. Also Silicone Injection Kits are usedfor supply the pressurized air to the silicone injector.

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Patient Population:

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Patients with various eye conditions.

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Contraindications:

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It should not be used outside of ophthalmic surgery.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training.Users must have operational experience.

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Factors That Can Affect the Success of the Operation:

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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The products are used with “EVA Anterior Machine NCBF” during surgery.

The purple apparatus (hose barb non valve insert) shown in the picture below is compatible with the “Dork” device if it is attached to the end of our product.

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If the white apparatus shown in the picture below is attached to the end of our product, it is compatible with the “Bausch Lomb” device.

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These parts are provided according to the request of the customers. These apparatuses are not sent to customers who want standard products.

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Accessory:

Aksesuar:

The product does not have any accessories.

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
Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

The packaging sample of our products is as seen in the photos below.

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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