Document No
Doküman No
TF08-MRN-IFU08-01

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

OPHTHALMIC SPATULA

OFTALMİK SPATULA

Company Name : MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve
Firma Adı
Authorized Representative : Cemil BAYTAŞ / General Manager
Firma Yetkilisi
Address : İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi
Adres
Phone : 0212 786 6253
Telefon
Web : www.merantip.com
Web

Product Models:

Ürün Modelleri:

Model and size information of the Ophthalmic Spatula products produced by us are shown in **TF08.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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There are two model of the product. They are “Lasik Flap Spatula” and “Spatula”.

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Lasik Flap Spatula Orucoglu;

Lastik Flap Spatula Oruçoğlu;

Lasik Flap Spatula is used to lift and reposition the cornea layer in eye surgery.

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Spatula;

Spatula;

Ophthalmic Spatulas are used for Cornea Scraper.

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Indications:

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Ophthalmic Spatulas are specially designed for various ophthalmic procedures. They are used for lifting or moving the cornea layer and for scraping on the cornea.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used outside of ophthalmic surgery.

It should not be used on infected eyes.

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Complications:

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After operation mild chemosis and redness can occur.

Intraoperative serous occurrence.

Iris laceration can occur. Lens capsule tear can occur.

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Revision Date
Revizyon Tarihi
20.11.2023

Side Effects:**Yan Etkiler:**

There are no known side effects of the product.

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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


- The membrane on the cornea layer is drawn with a spatula. Then the spatula enters between the membrane and the cornea and the membrane is completely removed from the cornea.

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Then, laser application is performed on the cornea layer.

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

- After this procedure, the membrane separated from the cornea is reassembled and the operation is completed by injecting various fluids into the eye area.

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Factors That Can Affect the Success of the Operation:

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Sterility/Sterilization Method:

Sterilite/Sterilizasyon Metodu:

The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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The products are used with no additional medical equipment.

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Accessory:**Aksesuar:**

The product does not have any accessories.

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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

Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

The packaging sample of our products is as seen in the photos below.

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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