Document No
Doküman No
TF08-MRN-IFU08-02

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

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Product Models:

Ürün Modelleri:

Model and size information of the Ophthalmic Manipulator products produced by us are shown in **TF08.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

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There are three model of the product. They are “IOL Manipulator Blunt Tip (Jaffe Maltzman)”, “IOL Manipulator Blunt Tip (Sinskey)” and “Nucleus Chopper”.

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IOL Manipulator Blunt Tip (Jaffe Maltzman):

IOL Manipulator Baskılama Ucu(Jaffe Maltzman)

It is used for placing Intraocular Lenses during cataract surgeries. It is not used for aspiration or irrigation, only for lens placement.

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It is used for placing Intraocular Lenses during cataract surgeries. It is not used for aspiration or irrigation, only for lens placement. Blunt tip avoids damage to intraocular lenses.

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Nucleus Chopper;

Nukleus Parçalama:

It is used for the manipulation of the cataracted eye lens / tissue in cataract surgery. The chopper tip is inclined back towards the nucleus to ensure the capture of even hard nuclei. May be used as a nucleus sustainer, the ball at the base of the chopper is smooth to protect the posterior capsule.

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Indications:

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Ophthalmic Manipulators are specially designed for various ophthalmic procedures. They are used for placing Intraocular Lenses during cataract surgeries. Also they are used for manipulation of the cataracted eye lens or tissue.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used outside of ophthalmic surgery. It should not be used on infected eyes.

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Revision Date
Revizyon Tarihi
20.11.2023

Complications:

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After operation mild chemosis and redness can occur.

Intraoperative serous occurrence.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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Usage Method:

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- During cataract surgery, with the help of various ophthalmic instruments, the cataracted tissue on the eye surface is broken. Then the Lens is put into the eye.

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


Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

- The lens placed in the eye is brought to the appropriate position with the Manipulators.

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- After the lens placement is completed, the manipulators are removed from the operation area and the cataract surgery procedure continues where it left off.

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Factors That Can Affect the Success of the Operation:

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Sterility/Sterilization Method:

Sterilite/Sterilizasyon Metodu:

The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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The products are used with no additional medical equipment.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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
Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

The packaging sample of our products is as seen in the photos below.

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