Document No
Doküman No
TF08-MRN-IFU08-03

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Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

MUSCLE HOOK

KAS KANCASI

Company Name
Firma Adı

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Authorized Representative
Firma Yetkilisi

: Cemil BAYTAŞ / General Manager

Address
Adres

: İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE

Phone
Telefon

: 0212 786 6253

Web
Web

: www.merantip.com

Product Models:

Ürün Modelleri:

Model and size information of the Muscle Hook products produced by us are shown in **TF08.TB.01**

Annex-1 Product Models Table with pictures.

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Intended Use:

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There are three model of the product. They are “Muscle Hook”, “Muscle Hook(Round)” and “Muscle Depressor(T)”.

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Muscle Hook:

Kas Kancası:

It is used to move the eye muscles up and down, right and left during the operation in eye surgery. In addition, Muscle Hooks in Eye Surgery are designed to elevate the muscle tissue for dissection and re-attachment in performing strabismus or other procedures with ocular muscles.

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Muscle Hook(Round):

Kas Kancası(Yuvarlak);

It is used to move the eyeball to the desired area during eye surgery or to ensure that the eyeball does not move during the operation.

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Muscle Depressor(T):

Kas Depresörü(T)

Muscle Depressor is used for controlling eye position. The product is used by the operator to depress the sclera of a patient's eye for either rotation or immobilization of the globe of the eye during examination.

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Indications:

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Muscle Hooks are generally used to move the eye muscles. In addition, it is used to provide or restrict the movement of the eyeball during the operation.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used outside of ophthalmic surgery.

It should not be used on infected eyes.

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Complications:

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After operation mild chemosis and redness can occur.

Intraoperative serous occurrence.

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Side Effects:

Yan Etkiler:

There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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

Usage Method:

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- It is used by the surgeon to control the movement of the eye muscles during the operation. After the necessary incisions are made on the eye surface, it is inserted under the muscle layer. It is directed by the surgeon to control the movements of the eye until the end of the operation.

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- After the operation is completed, these products are removed from the operation area.
Operasyon tamamlandıktan sonra bu ürünler operasyon alanından uzaklaştırılır.
- During the operation, the selection of the appropriate size and shape of the product should be made

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Factors That Can Affect the Success of the Operation:

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Sterility/Sterilization Method:

Sterilite/Sterilizasyon Metodu:

The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Shelf Life:

Raf Ömrü:

In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

Keep the product dry in any environment.

Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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The products are used with no additional medical equipment.

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Accessory:

Aksesuar:

The product does not have any accessories.

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Packaging, Labeling and Handling:

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Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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
Revision No
Revizyon No
04

Revision Date
Revizyon Tarihi
20.11.2023

The packaging sample of our products is as seen in the photos below.

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Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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