BACKFLUSH CANNULA
CANNULA BACKFLUSH

Company Name <i>Nome della società</i>	:	MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve Fatma Baytaş
Company Representative <i>Rappresentante dell'azienda</i>	:	Cemil BAYTAŞ / General Manager
Address <i>Indirizzo</i>	:	İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE
Phone <i>Telefono</i>	:	+90 212 786 6253
Web <i>Web</i>	:	www.merantip.com

Product Structure and General Features:

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The Backflush Cannula is designed for irrigation and aspiration to the operation area. Cannula tip is made of 304 Stainless Steel. The connection area at the back of the cannula is made of polycarbonate raw material. The products are available for single use and sterile. The cannulas to be used may differ depending on the operation area and the surgeon's choice. Models are available in various diameters and sizes in order to be suitable for every eye structure.

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Product Models:

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Model and size information of the Backflush Cannula products manufactured by us are shown in **TF03.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

Uso previsto:

Backflush Cannulas are used for irrigation and aspiration to the operation area. It also allows the irrigation and aspiration speed to be adjusted by the doctor and the control to be made by the doctor.

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Indications:

Indicazioni:

Backflush cannulas can be used in various procedures such as vacuum cleaning of sedimental blood from retinal surface, internal drainage of subretinal fluid, manipulation of flap, exchange of intraocular fluid for air or external liquids.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used outside of ophthalmic surgery.

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It should not be used on infected eyes.

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Complications:

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If the tip of the Cannula is not sharp, it may damage the tissue of the eye.

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After operation mild chemosis and redness can occur.

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Vitreous hemorrhage can occur.

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Intraoperative serous occurrence.

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Iatrogenic retinal breaks can occur.

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Side Effects:

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There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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Instruction for Use:

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Backflush Cannulas are restricted for use by surgeons trained in ophthalmic surgery.

- The Cannulas should be opened from the marked area without disruption of the sterility.

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- Cannula should be connected to the Equipment by a sterile person according to usage of the product.

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- After determining the area to be applied on, the cannula is inserted to exchange fluids or air actively or passively.

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- Dispose of used products in accordance with established hospital protocol for hazardous waste.

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Factors That Can Affect the Success Of The Operation:

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The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

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Keep the product dry in any environment.

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Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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The products are used with injectors during surgery.

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Accessory:

Accessorio:

The product does not have any accessories.

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Shelf Life:

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In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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


Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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The packaging sample of our products is as seen in the photos below.

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Symbols Prepared According to EN 15223-1 Standard and Their Meanings

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