LACRIMAL CANNULA
CANNULA LACRIMALE

Company Name
Nome della società : MMT MERAN MEDİKAL TEKNOLOJİ KOLL. ŞTİ. Ahmet Cemil Baytaş ve Fatma Baytaş

Company Representative
Rappresentante dell'azienda : Cemil BAYTAŞ / General Manager

Address
Indirizzo : İstanbul Trakya Serbest Bölgesi Ferhatpaşa S.B. Mah. Ali Rıza Efendi Cad. No:22/201 Çatalca/İSTANBUL/TÜRKİYE

Phone
Telefono : +90 212 786 6253

Web
Web : www.merantip.com

Product Structure and General Features:

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The Lacrimal Cannula is designed for the use in surgeries in the lacrimal duct and oculoplasty. Cannula tip is made of 304 Stainless Steel. The connection area at the back of the cannula is made of polycarbonate raw material. The products are available for single use and sterile. The cannulas to be used may differ depending on the operation area and the surgeon's choice. Models are available in various diameters and sizes in order to be suitable for every eye structure.

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Product Models:

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Model and size information of the Lacrimal Cannula products manufactured by us are shown in **TF03.TB.01 Annex-1 Product Models Table** with pictures.

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Intended Use:

Uso previsto:

Lacrimal Cannulas are disposable and used by surgeries in the lacrimal duct and oculoplasty. Lacrimal Cannulas also uses for the treatment of canaliculus system stenosis, dry eye, or nasolacrimal duct obstruction. These cannulas provide easy access directly to the irrigation and operation area.

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Indications:

Indicazioni:

Lacrimal Cannulas are used on eye capsules to tear and to mark eye capsule in cataract surgeries.

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Patient Population:

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Patients with various eye conditions, especially cataracts.

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Contraindications:

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It should not be used outside of ophthalmic surgery.

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It should not be used on infected eyes.

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Complications:

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If the tip of the Cannula is not sharp, it may damage the tissue of the eye.

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After operation mild chemosis and redness can occur.

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Vitreous hemorrhage can occur.

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Intraoperative serous occurrence.

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Side Effects:

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There are no known side effects of the product.

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Product Users:

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The product should be used by specialist physicians who have completed the necessary medical training. Users must have operational experience.

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Instruction for Use:

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- Lacrimal Cannulas are restricted for use by surgeons trained in ophthalmic surgery.

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- The Cannulas should be opened from the marked area without disruption of the sterility.

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- Cannula should be connected to the Equipment by a sterile person according to usage of the product.

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- Lacrimal Cannulas are disposable and used by surgeries in the lacrimal duct and oculoplasty. Lacrimal Cannulas also uses for the treatment of canaliculus system stenosis, dry eye, or nasolacrimal duct obstruction. These cannulas provide easy access directly to the irrigation and operation area.

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- Dispose of used products in accordance with established hospital protocol for hazardous waste.

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Factors That Can Affect the Success of The Operation:

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- The user's experience and experience in the operation,

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- Product selection in suitable type and size for the patient to be applied,

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- Eye structure of the patient,

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- Failure of the patient to have the necessary controls after the operation.

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The products are sterilized with Ethylene Oxide and presented to the user as sterile.

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The products are for single use only. There is absolutely no reusability.

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Product Storage and Shipping Conditions

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Protect the product against direct sunlight during shipment and storage.

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Keep the product dry in any environment.

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Store the products between 5 ° C and 35 ° C temperature and maximum 70% RH humidity.

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Medical Equipment Used with Products:

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The products are used with injectors during surgery.

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Accessory:

Accessorio:

The product does not have any accessories.

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Shelf Life:

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In line with the stability studies carried out in independent laboratories, the product shelf life was determined as 4 years.

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
Packaging material are manufactured according to EN ISO 11607-1 and EN 868-6 standard, tear, puncture and high resistance to liquids, which maintains sterility is the long packaging material. It sealed with heat cut in desired dimensions. Leak-proof barrier allows the product to remain sterile.

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The packaging sample of our products is as seen in the photos below.

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